


Product Datasheet

ISLR antibody (orb1132883)

Description	Rabbit polyclonal antibody to ISLR
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Raised against a recombinant human ISLR protein 111-275aa (BC022478).
Preservatives	PBS, pH 7.4 with 0.02% Sodium Azide
Form/Appearance	Liquid
Concentration	100ug/100ul
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Application notes	ELISA titer: 1:20,000-1:80,000 Western blot analysis: 1:500-1:1,000 Immunohistochemistry: 1:100-1:500 Optimal dilutions/concentrations should be determined by the end use.
Clonality	Polyclonal
Dilution Range	ELISA titer: 1:20,000-1:80,000 Western blot analysis: 1:500-1:1,000 Immunohistochemistry: 1:100-1:500 Optimal dilutions/concentrations should be determined by the end use.
Expiration Date	12 months from date of receipt.



WB analysis of fetal heart lysate using ...



IHC-P image of fetal heart tissue using ...