



Product Datasheet

Reelin antibody (orb11331)


Description	Reelin antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human Reelin. Please contact us for the exact immunogen sequence. The peptide is available as orb374744.
Target	Reelin
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
Uniprot ID	Q60841 , P78509 , P58751
NCBI	005045 , 21 , 31 , 005036
Entrez	5649
Dilution Range	IF/ICC: 1:100-800, IHC-P:1-100:800 (based on 0.5 mg/ml)
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of pig larg...



IHC-P image of 6pig large intestine tiss...



Immunohistochemical staining of mouse br...