

Product Datasheet

Ki67 antibody (orb1133242)

Description	Rabbit polyclonal antibody to Ki67
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant protein
Preservatives	PBS, pH 7.4 with 0.02% Sodium Azide
Form/Appearance	Liquid
Concentration	100ug/100ul
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Application notes	Optimal dilutions to be determined by researchers
Isotype	IgG
Clonality	Polyclonal
Dilution Range	Optimal dilutions to be determined by researchers
Expiration Date	12 months from date of receipt.