

Product Datasheet

ATP5B antibody (orb1133433)

Description

Rabbit polyclonal antibody to ATP5B

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications
Immunogen

Rabbit polyclonal ATP5B antibody was raised against a 16 amino acid peptide from near the intermediate of human ATP5B protein.

Preservatives

PBS, pH 7.4 with 0.02% Sodium Azide

Form/Appearance

Liquid

Concentration

100ug/100ul

Storage

This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Note

For research use only

Application notes

For research use only Peptide ELISA titer: 1:1,000-1:10,000 Western blot analysis: 1:500-1:1,000 Immunohistochemistry: 1:100-500 Optimal dilutions/concentrations should be determined by the end use

Isotype

IgG

Clonality

Polyclonal

Dilution Range

For research use only Peptide ELISA titer: 1:1,000-1:10,000 Western blot analysis: 1:500-1:1,000 Immunohistochemistry: 1:100-500 Optimal dilutions/concentrations should be determined by the end use

Expiration Date

12 months from date of receipt.