

Product Datasheet

DDX6 antibody (orb1133572)

Description	Rabbit polyclonal antibody to DDX6
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthetic peptide conjugated to KLH
Preservatives	PBS, pH 7.4 with 0.02% Sodium Azide
Form/Appearance	Liquid
Concentration	100ug/100ul
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Application notes	For research use only Peptide ELISA: 1:40,000-80,000 Western blot analysis: 1:200-500 Optimal dilutions/concentrations should be determined by the end use
Isotype	Rabbit IgG
Clonality	Polyclonal
Dilution Range	For research use only Peptide ELISA: 1:40,000-80,000 Western blot analysis: 1:200-500 Optimal dilutions/concentrations should be determined by the end use
Expiration Date	12 months from date of receipt.