

## Product Datasheet

### APE1 antibody (orb1133633)

<b>Description</b>	Rabbit polyclonal antibody to APE1
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, WB
<b>Immunogen</b>	Synthetic peptide conjugated to KLH
<b>Preservatives</b>	PBS, pH 7.4 with 0.02% Sodium Azide
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	100ug/100ul
<b>Storage</b>	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
<b>Note</b>	For research use only
<b>Application notes</b>	Recommended dilutions Peptide ELISA: 1:80,000 WB: 1/200 - 1/1000 Not yet tested in other applications. Optimal dilutions/concentrations should be determined by the end user
<b>Isotype</b>	Rabbit IgG
<b>Clonality</b>	Polyclonal
<b>Dilution Range</b>	Recommended dilutions Peptide ELISA: 1:80,000 WB: 1/200 - 1/1000 Not yet tested in other applications. Optimal dilutions/concentrations should be determined by the end user
<b>Expiration Date</b>	12 months from date of receipt.