

Product Datasheet

VDAC1 antibody (orb1133660)

Description

Rabbit polyclonal antibody to VDAC1

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, WB

Applications
Immunogen

Synthetic peptide conjugated to KLH

Preservatives

PBS, pH 7.4 with 0.02% Sodium Azide

Form/Appearance

Liquid

Concentration

100ug/100ul

Storage

This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Note

For research use only

Application notes

Recommended dilutions Peptide ELISA: 1:20,000 WB: 1/100 - 1/500 Not yet tested in other applications. Optimal dilutions/concentrations should be determined by the end user

Isotype

Rabbit IgG

Clonality

Polyclonal

Dilution Range

Recommended dilutions Peptide ELISA: 1:20,000 WB: 1/100 - 1/500 Not yet tested in other applications. Optimal dilutions/concentrations should be determined by the end user

Expiration Date

12 months from date of receipt.