

Product Datasheet

SARS orf1ab polyprotein antibody (orb11349)

Description	SARS orf1ab polyprotein antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from SARS putative orflab polyprotein (3501-3600/7073 aa)
Target	SARS putative orflab polyprotein
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	790 kDa
Uniprot ID	P59594
Expiration Date	12 months from date of receipt.