



Product Datasheet

SATB1 antibody (orb11351)

Description	SATB1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	KLH conjugated synthetic peptide derived from mouse SATB1 (671-764/764 aa)
Target	SATB1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	86 kDa
Uniprot ID	Q01826
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...



IHC-P of human tumor tissue (SATB1 antib...