

## Product Datasheet

### CXCL14 antibody (orb11354)

## Description

CXCL14 antibody

### Species/Host

Rabbit

### Reactivity

Human, Mouse, Rat

### Conjugation

Unconjugated

### Tested

ELISA, IHC-P, WB

### Applications

### Immunogen

KLH conjugated synthetic peptide derived from human CXCL14. Please contact us for the exact immunogen sequence. The peptide is available as orb374746.

### Target

CXCL14

### Form/Appearance

10 mM PBS, 0.02% sodium azide

### Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

### Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

### Note

For research use only

### Isotype

IgG

### Clonality

Polyclonal

### Purity

Polyclonal antibodies are purified by peptide affinity chromatography

### MW

13 kDa

### Uniprot ID

[O95715](#)

### Entrez



[9547](#)

### Dilution Range


WB: 1:100-2000, IHC-P: 1:50-300 (based on 0.5 mg/ml)

### Expiration Date


12 months from date of receipt.



Immunohistochemical staining of rat brain...



IHC-P staining of rat brain tissue using...



IHC-P staining of rat pancreas tissue us...