

## Product Datasheet

**Tetanus Toxin fragment C antibody [Out of stock]  
(orb1136274)**

## Description

Rabbit polyclonal antibody to

### Species/Host

Rabbit

### Reactivity

Other

### Conjugation

Unconjugated

### Tested

ELISA, IHC-P, WB

### Applications

### Immunogen

This antibody was prepared from whole rabbit serum produced by repeated immunizations with a polyhistidine-tagged TTFC fusion protein expressed in E. coli corresponding to full length TTFC. The polyhistidine-tag was placed at the amino terminus of the protein. The recombinant protein was purified from soluble bacterial extracts by metal-affinity and ion-exchange chromatography.

### Preservatives

PBS, pH 7.4 with 0.05% sodium azide.

### Form/Appearance

Liquid

### Concentration

1ug/ul

### Storage

This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

### Note

For research use only

### Isotype

IgG

### Clonality

Polyclonal

### Purity

Purified

### Dilution Range

This antibody has been tested for use in ELISA and western blot. Specific conditions for reactivity should be optimized by the end user. For western blot experiments, one can expect a band of  $\approx 52$  kDa in size corresponding to full-length

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com) | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive<br>Research Triangle Park<br>Durham, North Carolina<br>27709, United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com) | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

ELISA, 0.1 ml of the undiluted antiserum (1:4000) was able to detect 6