

Product Datasheet

Choline Acetyltransferase (ChAT) antibody [Out of stock] (orb1137621)

Description	Rabbit polyclonal antibody to Choline Acetyltransferase (ChAT)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	22 amino acid peptide from porcine ChAT (GLF SSY RLP GHT QDT LVA QKSS) coupled to KLH
Preservatives	PBS, pH 7.4 with 0.05% sodium azide.
Form/Appearance	Liquid
Concentration	1ug/ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purified
Expiration Date	12 months from date of receipt.