

Product Datasheet

**Glutamine Synthetase antibody [Out of stock]
(orb1137719)**

Description

Goat polyclonal antibody against

	Biotin
Species/Host	Goat
Reactivity	Other
Conjugation	Biotin
Tested Applications	WB
Immunogen	Glutamine Synthetase [Microbial]
Preservatives	PBS, pH 7.4 with 0.05% sodium azide.
Form/Appearance	Liquid
Concentration	1ug/ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Ion exchange chromatography
Dilution Range	Suitable for immunoblotting, ELISA, immunohistochemistry, immunomicroscopy as well as other antibody based assays using streptavidin or avidin conjugates requiring lot-to-lot consistency.
Expiration Date	12 months from date of receipt.