

Product Datasheet

Sortilin (Neurotensin receptor 3) antibody (orb1138063)

Description

Rabbit polyclonal antibody to Sortilin

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Peptide DKDTTRRIHVSTDQD(C) corresponding to residues 320-335 of human Sortilin (Accession Q99523).
Preservatives	PBS, pH 7.4 with 0.05% sodium azide.
Form/Appearance	Liquid
Concentration	1ug/ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Application notes	For research purposes only, not for human use.
Isotype	IgG
Clonality	Polyclonal
Dilution Range	For research purposes only, not for human use.
Expiration Date	12 months from date of receipt.