

Product Datasheet

Nogo Receptor (Nogo-R, Reticulon receptor 4, NgR1, Nogo-66 Receptor) antibody (orb1138070)

Description

Rabbit polyclonal antibody to Nogo

Species/Host

Rabbit

Reactivity

Mouse, Rat

Conjugation

Unconjugated

Tested

ICC, IHC, WB

Applications
Immunogen

Peptide (C)GDSPPGNGSGPRHIND corresponding to amino acid residues 366-381 of the human Nogo receptor (Accession Q9BZR6).

Preservatives

PBS, pH 7.4 with 0.05% sodium azide.

Form/Appearance

Liquid

Concentration

1ug/ul

Storage

This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Note

For research use only

Application notes

For research purposes only, not for human use.

Isotype

IgG

Clonality

Polyclonal

Dilution Range

For research purposes only, not for human use.

Expiration Date

12 months from date of receipt.