

Product Datasheet

Cannabinoid Receptor 1 (extracellular) antibody
[Out of stock] (orb1138105)

Description

Rabbit polyclonal antibody to

(extracellular)

Species/Host

Rabbit

Reactivity

Mouse, Rat

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications
Immunogen

Peptide NKSLSSFKENEENIQ
corresponding to amino acid
residues 84-99 of rat CB1 receptor
(Accession P20272).

Preservatives

PBS, pH 7.4 with 0.05% sodium
azide.

Form/Appearance

Liquid

Concentration

1ug/ul

Storage

This product is stable for several
weeks at 4° C as an undiluted
liquid. Dilute only prior to
immediate use. For extended
storage, aliquot contents and
freeze at -20° C or below. Avoid
cycles of freezing and thawing.
Expiration date is one (1) year
from date of receipt.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Affinity purified on immobilized
antigen.

Dilution Range

For research purposes only, not for
human use.

Expiration Date

12 months from date of receipt.