

Product Datasheet

ADAM-15 (A Disintegrin And Metalloproteinase-15, Metargidin, MDC15), Cytoplasmic domain antibody [Out of stock] (orb1138274)

Description

Rabbit polyclonal antibody to ADAM-15

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

WB

Applications
Immunogen

Synthetic peptide based on the cytoplasmic domain of human ADAM15.

Preservatives

PBS, pH 7.4 with 0.05% sodium azide.

Form/Appearance

Liquid

Concentration

1ug/ul

Storage

This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Immunogen Affinity purified

Dilution Range

Recommended dilutions WB:
Recommended starting dilution of 1/1000 (when using colorimetric substrates such as BCIP/NBT) and 1/5000 (for chemiluminescent substrates). Detects a band of approximately 90 kDa (predicted molecular weight: 88 kDa). Note: The antibody recognizes bands of 110 kDa and 90 kDa, the zymogen and furin processed forms of ADAM15. A series of cleaved products recover at 55 kDa and 40 kDa can also be seen in reduced Western blots of cell lysates. Dilution optimised using Chromogenic detection. Not yet tested in other