

Product Datasheet

**Estrogen Receptor- beta antibody [Out of stock]
(orb1138695)**

Description

Mouse monoclonal antibody to

Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Ni-NTA purified recombinant human ER-alpha (410-592aa) expressed in E. Coli strain BL21 (DE3)
Preservatives	PBS, pH 7.4 with 0.05% sodium azide.
Form/Appearance	Liquid
Concentration	1ug/ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	purified by protein A affinity chromatography.
Dilution Range	Western Blot: Dilution 1: 100- 1: 500; IHC (P): Dilution 1: 100- 1: 500; IHC (F): Dilution 1: 100- 1: 500; Determining optimal working dilutions by titration test.
Expiration Date	12 months from date of receipt.