

Product Datasheet

GFP(Green fluorescent protein) antibody (orb1138727)

Description	Rabbit polyclonal antibody to GFP(Green fluorescent protein)
Species/Host	Rabbit
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Polyclonal antibody produced in rabbits immunizing with a synthetic peptide corresponding to residues 70-83 of Aequorea victoria GFP (Green fluorescent protein)
Preservatives	PBS, pH 7.4 with 0.05% sodium azide.
Form/Appearance	Liquid
Concentration	1ug/ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.