

Product Datasheet

CHT1(High-affinity choline transporter 1) antibody
[Out of stock] (orb1138999)

Description

Rabbit polyclonal antibody to

transporter 1)

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications
Immunogen

Polyclonal antibody produced in rabbits immunizing with a synthetic peptide corresponding to C-terminal residues of human choline transporter

Preservatives

PBS, pH 7.4 with 0.05% sodium azide.

Form/Appearance

Liquid

Concentration

1ug/ul

Storage

This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Purified by antigen-specific affinity chromatography.

Expiration Date

12 months from date of receipt.