

Product Datasheet

Estrogen Receptor- beta (Phospho-Ser118) antibody (orb1140150)

Description

Rabbit polyclonal antibody to

Ser118)

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

IHC, WB

Applications
Immunogen

The antiserum was produced against synthesized phosphopeptide derived from human Estrogen Receptor- α around the phosphorylation site of serine 118 (Q-L-SP-P-F).

Preservatives

phosphate buffered saline (without Mg^{2+} and Ca^{2+}), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Form/Appearance

Liquid

Concentration

1 μ g/ μ l

Storage

This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Note

For research use only

Application notes

WB: 1:500~1:1000 IHC: 1:50~1:100

Clonality

Polyclonal

Dilution Range

WB: 1:500~1:1000 IHC: 1:50~1:100

Expiration Date

12 months from date of receipt.