

Product Datasheet

Recombinant Ki-67 antibody (orb1140240)

Description	Rabbit monoclonal antibody to Recombinant Ki-67
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Protein
Form/Appearance	Tris Buffer, pH 7.3 - 7.7, with 1% rAlbumin and 0.1% Sodium Azide
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Application notes	1:50-1:200
Isotype	IgG, kappa
Clonality	Monoclonal
Dilution Range	IHC: 1:1000-1:3000; ELISA: 1:1000-1:2000 (Optimal dilutions or concentrations should be determined by the end user)
Expiration Date	12 months from date of receipt.