

Product Datasheet

CFP Tag antibody (orb1140373)

Description	Mouse monoclonal antibody to CFP Tag
Species/Host	Mouse
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	ELISA, IP, WB
Immunogen	GFP from the jellyfish Aequorea victoria N-Terminal peptide-KLH conjugates (first 20 aa)
Form/Appearance	Liquid
Concentration	100ug/100ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Application notes	WB 1:1000. For other applications, the researcher should determine the optimal working dilution.
Isotype	IgG
Clonality	Monoclonal
Dilution Range	WB 1:1000. For other applications, the researcher should determine the optimal working dilution.
Expiration Date	12 months from date of receipt.