

Product Datasheet

alpha-TubulinLoading control antibody (orb1140402)

Description

Mouse monoclonal antibody to alpha-

Species/Host	Mouse
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	Recombinant α -tubulin protein from bacteria.
Form/Appearance	Liquid
Concentration	100ug/100ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Application notes	Western blot at dilution of 1:1000 and Immunofluorescence staining at 1:100 dilutions. It is very useful as a loading control for quantitative western blotting. The antibody recognizes a 50kDa band in a Western Blot analysis.
Isotype	IgG
Clonality	Monoclonal
Dilution Range	Western blot at dilution of 1:1000 and Immunofluorescence staining at 1:100 dilutions. It is very useful as a loading control for quantitative western blotting. The antibody recognizes a 50kDa band in a Western Blot analysis.
Expiration Date	12 months from date of receipt.