

Product Datasheet

alpha Spectrin antibody (orb11408)

Description	alpha Spectrin antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	IF
Immunogen	KLH conjugated synthetic peptide derived from human Spectrin alpha chain, non-erythrocytic 1 (30-110/2472 aa)
Target	Spectrin alpha chain, non-erythrocytic 1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	272 kDa
Uniprot ID	Q13813
Dilution Range	IF=1:100-500, IHC-P=1:100-500, IHC-F=1:100-500, ICC=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.