



Product Datasheet

SEC14L2 antibody (orb11410)

Description	Rabbit polyclonal antibody to SEC14L2
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human SEC14L2(NP_001191133.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	46kDa
Uniprot ID	O76054
Dilution Range	WB: 1:1000-1:2000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U-2 OS cells...



Western blot analysis of extracts of various tissues...