

## Product Datasheet

# Rabbit anti-COVID-19/SARS-CoV-2 S1-S2 Polyclonal Antibody (orb1141353)

|                            |                                                           |
|----------------------------|-----------------------------------------------------------|
| <b>Description</b>         | Rabbit anti-COVID-19/SARS-CoV-2 S1-S2 Polyclonal Antibody |
| <b>Species/Host</b>        | Rabbit                                                    |
| <b>Endotoxins</b>          | < 0.1EU/ug determined by LAL method.                      |
| <b>Conjugation</b>         | Unconjugated                                              |
| <b>Tested Applications</b> | ELISA                                                     |
| <b>Note</b>                | For research use only                                     |
| <b>Clonality</b>           | Polyclonal                                                |
| <b>Expiration Date</b>     | 12 months from date of receipt.                           |