

Product Datasheet

Recombinant Human anti-COVID-19/SARS-CoV-2 S1 Monoclonal Antibody (orb1141359)

Description	Recombinant Human anti-COVID-19/SARS-CoV-2 S1 Monoclonal Antibody
Species/Host	Human
Endotoxins	< 0.1EU/ug determined by LAL method.
Conjugation	Unconjugated
Tested Applications	Control, ELISA, GICA, NeA
Note	For research use only
Isotype	IgM
Clonality	Recombinant
Expiration Date	12 months from date of receipt.