

Product Datasheet

Recombinant Human anti-COVID-19/SARS-CoV-2 Nucleocapsid Protein Monoclonal Antibody (orb1141367)

Description	Recombinant Human anti-COVID-19/SARS-CoV-2 Nucleocapsid Protein Monoclonal Antibody
Species/Host	Human
Endotoxins	< 0.1EU/ug determined by LAL method.
Conjugation	Unconjugated
Tested Applications	ELISA, GICA, WB
Note	For research use only
Isotype	IgG1
Clonality	Recombinant
Expiration Date	12 months from date of receipt.