



Product Datasheet

SRF antibody (orb11414)


Description	SRF antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human SRF (101-250/508 aa)
Target	SRF
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	56 kDa
Uniprot ID	P31948
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, ICC=1:50, Flow-Cyt=1ug/Test
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse intestine...



Western blot analysis of mouse lysates u...



Western blot analysis of Mouse Muscle ly...