

Product Datasheet

Recombinant SARS-CoV-2/COVID-19 PL-Pro (orb1141418)

Description	Recombinant SARS-CoV-2/COVID-19 PL-Pro
Species/Host	E. coli
Note	For research use only
Purity	≥ 90% as determined by SDS-PAGE
Source	Recombinant SARS-CoV-2 PL-Pro was expressed from E.coli expression system and the target gene encoding Glu1564-Ile1877 is expressed with N- His Tag
Expiration Date	6 months from date of receipt.