

Product Datasheet

Mouse Interferon gamma, IFN-gamma ELISA Kit (orb1141470)

Description	Mouse Interferon γ ,IFN- γ ELISA Kit
Reactivity	Mouse
Range	0.156 ng/mL-10 ng/mL
Target	IFNG
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P01580
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.