

Product Datasheet

Rat fibroblast growth factor 23 (FGF23) ELISA Kit (orb1141474)

Description	Rat fibroblast growth factor 23 (FGF23) ELISA Kit
Reactivity	Rat
Range	15.6 ng/mL-1000 ng/mL
Target	FGF23
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	Q8VI82
Sensitivity	3.9 ng/mL
Expiration Date	Please enquire.