

Product Datasheet

IGFBP-1 antibody (orb1141512)

Description	Mouse monoclonal antibody to IGFBP-1
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	CLIA, EIA, LF
Target	IGFBP-1
Preservatives	0.05% Proclin 300. PBS, pH 7.4.
Form/Appearance	Supplied as a 0.22 µm filtered solution of PBS, pH 7.4
Storage	Short term (up to 7 days) store at 2-8°C. Long term, store at -20°C. Avoid multiple freeze/thaw cycles.
Note	For research use only
Application notes	Method: Sandwich method. Use: Coating
Isotype	IgG
Clonality	Monoclonal
Purity	> 95% as determined by SDS-PAGE.
Uniprot ID	P08833
Expiration Date	12 months from date of receipt.