

Product Datasheet

WASL antibody (orb1141629)

Description	Mouse monoclonal antibody to WASL
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human WASL (NP_003932) produced in HEK293T cell.
Preservatives	PBS (PH 7.3) containing 1% rAlbumin, 50% glycerol and 0.02% sodium azide.
Concentration	1mg/ml
Storage	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Note	For research use only
Clonality	Monoclonal
MW	54.6 kDa(Predicted)
Uniprot ID	O00401
Dilution Range	WB 1:500~2000
Expiration Date	12 months from date of receipt.