

Product Datasheet

UQCRC1 antibody (orb1141675)

Description	Mouse monoclonal antibody to UQCRC1
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 46-305 of human UQCRC1 (NP_003356) produced in E.coli.
Preservatives	PBS (PH 7.3) containing 1% rAlbumin, 50% glycerol and 0.02% sodium azide.
Concentration	1mg/ml
Storage	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Note	For research use only
Clonality	Monoclonal
MW	52.5 kDa(Predicted)
Uniprot ID	P31930
Dilution Range	WB 1:500~1000, IHC 1:150, FLOW 1:100
Expiration Date	12 months from date of receipt.