

## Product Datasheet

### SORL1 antibody (orb1142085)

|                            |                                                                                                                        |
|----------------------------|------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to SORL1                                                                                     |
| <b>Species/Host</b>        | Mouse                                                                                                                  |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                                                      |
| <b>Conjugation</b>         | Unconjugated                                                                                                           |
| <b>Tested Applications</b> | IHC, WB                                                                                                                |
| <b>Immunogen</b>           | Human recombinant protein fragment corresponding to amino acids 573-753 of human SORL1 (NP_003096) produced in E.coli. |
| <b>Preservatives</b>       | PBS (PH 7.3) containing 1% rAlbumin, 50% glycerol and 0.02% sodium azide.                                              |
| <b>Concentration</b>       | 1mg/ml                                                                                                                 |
| <b>Storage</b>             | PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.                                                   |
| <b>Note</b>                | For research use only                                                                                                  |
| <b>Clonality</b>           | Monoclonal                                                                                                             |
| <b>MW</b>                  | 248.4 kDa(Predicted)                                                                                                   |
| <b>Uniprot ID</b>          | <b>Q92673</b>                                                                                                          |
| <b>Dilution Range</b>      | WB 1:500~2000, IHC 1:2000                                                                                              |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                        |