

Product Datasheet

SH3PXD2A antibody (orb1142163)

Description	Mouse monoclonal antibody to SH3PXD2A
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 299-395 of human SH3PXD2A (NP_055446) produced in E.coli.
Preservatives	PBS (PH 7.3) containing 1% rAlbumin, 50% glycerol and 0.02% sodium azide.
Concentration	1mg/ml
Storage	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Note	For research use only
Clonality	Monoclonal
MW	121.7 kDa(Predicted)
Uniprot ID	B3KPL1
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.