

Product Datasheet

SETD2 antibody (orb1142181)

Description	Mouse monoclonal antibody to SETD2
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1787-2144 of human SETD2 (NP_054878) produced in E.coli.
Preservatives	PBS (PH 7.3) containing 1% rAlbumin, 50% glycerol and 0.02% sodium azide.
Concentration	1mg/ml
Storage	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Note	For research use only
Clonality	Monoclonal
MW	287 kDa(Predicted)
Uniprot ID	Q9BYW2
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.