



Product Datasheet

SSTR2 antibody (orb11422)

Description	SSTR2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	KLH conjugated synthetic peptide derived from human SSTR2 (211-320/369 aa)
Target	Somatostatin Receptor 2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	41 kDa
Uniprot ID	P30874
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...



IHC-P of mice lung tissue (SSTR2 antibody)