

Product Datasheet

RPA2 antibody (orb1142323)

Description	Mouse monoclonal antibody to RPA2
Species/Host	Mouse
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, IP, WB
Immunogen	Full length human recombinant protein of human RPA2 (NP_002937) produced in HEK293T cell.
Preservatives	PBS (pH 7.3) containing 1% rAlbumin, 50% glycerol and 0.02% sodium azide.
Concentration	1mg/ml
Storage	PBS (pH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Note	For research use only
Clonality	Monoclonal
MW	29.1 kDa(Predicted)
Uniprot ID	B4DUL2
Dilution Range	WB 1:1000~2000, IHC 1:50, IF 1:100, FLOW 1:100, IP 2ug/500ul
Expiration Date	12 months from date of receipt.