

Product Datasheet

RFC3 antibody (orb1142395)

Description	Mouse monoclonal antibody to RFC3
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human RFC3 (NP_002906) produced in E.coli.
Preservatives	PBS (PH 7.3) containing 1% rAlbumin, 50% glycerol and 0.02% sodium azide.
Concentration	1mg/ml
Storage	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Note	For research use only
Clonality	Monoclonal
MW	40.4 kDa(Predicted)
Uniprot ID	A0A024RDQ8
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.