

Product Datasheet

RFC2 antibody (orb1142399)

Description	Mouse monoclonal antibody to RFC2
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of huamn RFC2(NP_002905) produced in E.coli.
Preservatives	PBS (PH 7.3) containing 1% rAlbumin, 50% glycerol and 0.02% sodium azide.
Concentration	1mg/ml
Storage	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Note	For research use only
Clonality	Monoclonal
MW	35.1 kDa(Predicted)
Uniprot ID	P35250
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.