

Product Datasheet

LRP6 antibody (Biotin) (orb114240)

Description	Rabbit polyclonal to LRP6 (Biotin).
Species/Host	Rabbit
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	KLH conjugated synthetic peptide derived from human LRP6
Target	LRP6
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Application notes	orb114240 is the Biotin conjugated form of orb6329
Isotype	IgG
Clonality	Polyclonal
MW	175kDa
Uniprot ID	O75581
Expiration Date	12 months from date of receipt.