

## Product Datasheet

### RET antibody (orb1142402)

|                            |                                                                                                                      |
|----------------------------|----------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to RET                                                                                     |
| <b>Species/Host</b>        | Mouse                                                                                                                |
| <b>Reactivity</b>          | Human, Rat                                                                                                           |
| <b>Conjugation</b>         | Unconjugated                                                                                                         |
| <b>Tested Applications</b> | IHC, WB                                                                                                              |
| <b>Immunogen</b>           | Human recombinant protein fragment corresponding to amino acids 713-1017 of human RET(NP_066124) produced in E.coli. |
| <b>Preservatives</b>       | PBS (PH 7.3) containing 1% rAlbumin, 50% glycerol and 0.02% sodium azide.                                            |
| <b>Concentration</b>       | 1mg/ml                                                                                                               |
| <b>Storage</b>             | PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.                                                 |
| <b>Note</b>                | For research use only                                                                                                |
| <b>Clonality</b>           | Monoclonal                                                                                                           |
| <b>MW</b>                  | 124.1 kDa(Predicted)                                                                                                 |
| <b>Uniprot ID</b>          | <b>P07949</b>                                                                                                        |
| <b>Dilution Range</b>      | IHC 1:150                                                                                                            |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                      |