

Product Datasheet

PIK3CG antibody (orb1142793)

Description	Mouse monoclonal antibody to PIK3CG
Species/Host	Mouse
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human PIK3CG(NP_002640) produced in HEK293T cell.
Preservatives	PBS (PH 7.3) containing 1% rAlbumin, 50% glycerol and 0.02% sodium azide.
Concentration	1mg/ml
Storage	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Note	For research use only
Clonality	Monoclonal
MW	126.3 kDa(Predicted)
Uniprot ID	P48736
Dilution Range	WB 1:200~2000, IHC 1:150, IF 1:100
Expiration Date	12 months from date of receipt.