

Product Datasheet

PIK3C2A antibody (orb1142807)

Description	Mouse monoclonal antibody to PIK3C2A
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 230-560 of human PIK3C2A (NP_002636) produced in E.coli.
Preservatives	PBS (PH 7.3) containing 1% rAlbumin, 50% glycerol and 0.02% sodium azide.
Concentration	1mg/ml
Storage	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Note	For research use only
Clonality	Monoclonal
MW	190.5 kDa(Predicted)
Uniprot ID	O00443
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.