

Product Datasheet

PFDN6 (PFD6) antibody (orb1142834)

Description	Mouse monoclonal antibody to PFDN6 (PFD6)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 2-129 of human PFDN6(NP_055075) produced in E.coli.
Preservatives	PBS (PH 7.3) containing 1% rAlbumin, 50% glycerol and 0.02% sodium azide.
Concentration	1mg/ml
Storage	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Note	For research use only
Clonality	Monoclonal
MW	14.4 kDa(Predicted)
Uniprot ID	Q5STK2
Dilution Range	WB 1:2000, IHC 1:150, FLOW 1:100
Expiration Date	12 months from date of receipt.