

Product Datasheet

PCBP1 antibody (orb1142918)

Description	Mouse monoclonal antibody to PCBP1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-225 of human PCBP1(NP_006187) produced in E.coli.
Preservatives	PBS (PH 7.3) containing 1% rAlbumin, 50% glycerol and 0.02% sodium azide.
Concentration	1mg/ml
Storage	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Note	For research use only
Clonality	Monoclonal
MW	37.3 kDa(Predicted)
Uniprot ID	Q53SS8
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.