

Product Datasheet

NTRK1 antibody (orb1143025)

Description	Mouse monoclonal antibody to NTRK1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 440-796 of human NTRK1(NP_002520) produced in E.coli.
Preservatives	PBS (PH 7.3) containing 1% rAlbumin, 50% glycerol and 0.02% sodium azide.
Concentration	1mg/ml
Storage	Store at -20C° for 2 years
Note	For research use only
Clonality	Monoclonal
MW	87.3 kDa(Predicted)
Uniprot ID	P04629
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.